


Anti-Ndufb7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc - P17568

Item Number: CM0014215



Introduction

Rabbit polyclonal antibody against human NDUFB7 (NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 7) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat mitochondrial complex I subunit B18. Supports mitochondrial disease research and IVD assay development.

[Learn More](#)

Parameter	Specification
Product Name	Anti-NDUFB7 Rabbit Polyclonal Antibody
Remarks/Alias	B18, CI-B18, MC1DN39, NDUFB7
Species	Human
Gene ID (Human)	4713
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-137 of human NDUFB7 (NP_004137.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17568
Protein Weight	16kDa
Shipping	Ice bag