


Anti-Ndufb8 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95169

Item Number: CM0024201



Introduction

High-quality rabbit polyclonal antibody against human NDUFB8, validated for WB, IHC-P, IF/ICC, ELISA. Specific for human, mouse, rat mitochondrial complex I accessory subunit. Ideal for mitochondrial disease research.

[Learn More](#)

Parameter	Description
Product Name	NDUFB8 Rabbit pAb
Remarks/Alias	ASH1, CI-ASH1, MC1DN32, NDUFB8
Species	Human
GeneID (Human)	4714
GeneID	4714
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-186 of human NDUFB8 (NP_004995.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95169
Protein Weight	22kDa
Shipping	Ice bag