

Anti-Ndufb8 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - O95169

Item Number: CM0002115



Introduction

Recombinant rabbit monoclonal antibody against human NDUFB8, a mitochondrial Complex I accessory subunit. Validated in WB, IF/ICC, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for IVD raw material supply and mitochondrial dysfunction research.

[Learn More](#)

Parameter	Value
Product Name	NDUFB8 Rabbit mAb
Remarks/Alias	ASHI; CI-ASHI; MC1DN32; NDUFB8
Species	Human
Gene ID (Human)	4714
Gene ID	4714
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human NDUFB8 (O95169).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95169
Protein Weight	22kDa
Shipping	Ice bag