

Anti-Ndufb9 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6M9

Item Number: CM0002250



Introduction

Rabbit monoclonal antibody against human NDUF9, a mitochondrial Complex I subunit. Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUF9 Rabbit mAb
Remarks/Alias	B22; LYRM3; CI-B22; UQOR22; NDUF9
Species	Human
GeneID (Human)	4715
GeneID	4715
Immunogen	Recombinant protein corresponding to amino acids 80-170 of human NDUF9 (Q9Y6M9)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6M9
Protein Weight	22kDa
Shipping	Ice bag