


Anti-Ndufc2 Polyclonal Antibody For Wb, Elisa - 095298

Item Number: CM0029564



Introduction

Rabbit polyclonal antibody targeting NDUF2 (Complex I subunit B14.5b) for WB and ELISA applications. Cross-reacts with human, mouse, rat. Ideal for mitochondrial respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUFC2 Rabbit pAb
Remarks/Alias	HLC-1; B14.5b; MC1DN36; NADHDH2; CI-B14.5b; NDUFC2
Species	Human
Gene ID (Human)	4718
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human NDUFC2 (NP_004540.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95298
Protein Weight	14kDa
Shipping	Ice bag