


Anti-Ndufc2 Rabbit Monoclonal Antibody For Wb, Elisa - O95298

Item Number: CM0002440



Introduction

Rabbit monoclonal antibody against human NDUFC2 (CI-B14.5b), validated for WB and ELISA. Targets the N-terminal region (aa 1-119) of this 14 kDa mitochondrial complex I accessory subunit, supporting research on electron transport chain assembly and mitochondrial diseases.

[Learn More](#)

Parameter	Value
Product Name	NDUFC2 Rabbit mAb
Remarks/Alias	HLC-1; B14.5b; MC1DN36; NADHDH2; CI-B14.5b; NDUFC2
Species	Human
GeneID (Human)	4718
GeneID	4718
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1-119 of human NDUFC2 (NP_004540.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95298
Protein Weight	14kDa
Shipping	Ice bag