


Anti-Ndufs1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P28331

Item Number: CM0007465



Introduction

Rabbit monoclonal antibody targeting human NDUFS1, suited for WB, IHC-P, and ELISA. Detects endogenous NDUFS1 in human, mouse, and rat samples. Designed for mitochondrial complex I studies and respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUFS1 Rabbit mAb
Remarks/Alias	CI-75k; MC1DN5; CI-75Kd; PRO1304; NDUFS1
Species	Human
GeneID (human)	4719
GeneID	4719
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 80-290 of human NDUFS1 (NP_004997.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28331
Protein Weight	79kDa
Shipping	Ice bag