

Anti-Ndufs1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P28331

Item Number: CM0029683



Introduction

NDUF51 Rabbit pAb validated for WB, IHC-P, ELISA. Targets human NDUF51 (CI-75kD), a core Complex I subunit essential for mitochondrial electron transport. Cross-reacts with mouse, rat. Ideal for mitochondrial disease research.

[Learn More](#)

Product Name	NDUF51 Rabbit pAb
Remarks/Alias	CI-75k; MC1DN5; CI-75Kd; PRO1304; NDUF51
Species	Human
GeneID (Human)	4719
GeneID	4719
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-290 of human NDUF51 (NP_004997.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28331
Protein Weight	79kDa
Shipping	Ice bag