

Anti-Ndufs1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28331

Item Number: CM0023940



Introduction

NDUFS1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Targets human NADH-ubiquinone oxidoreductase 75 kDa subunit (Complex I-75kD). Cross-reacts with mouse, rat. Ideal for mitochondrial complex I research and diagnostics.

[Learn More](#)

Parameter	Value
Product Name	NDUFS1 Rabbit pAb
Remarks/Alias	CI-75k, MC1DN5, CI-75Kd, PRO1304, NDUFS1
Species	Human
GeneID (Human)	4719
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-290 of human NDUFS1 (NP_004997.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28331
Protein Weight	79kDa
Shipping	Ice bag