


Anti-Ndufs2 Monoclonal Antibody For Wb And Elisa - 075306

Item Number: CM0003501



Introduction

Rabbit monoclonal antibody targeting human NDUFS2 (CI-49), suitable for WB and ELISA. Cross-reacts with mouse and rat. Validated for research on mitochondrial Complex I and oxidative phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	NDUFS2 Rabbit mAb
Remarks/Alias	CI-49; MC1DN6; NDUFS2
Species	Human
Gene ID	4720
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 300-434 of human NDUFS2 (NP_004541.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75306
Protein Weight	53kDa
Shipping	Ice bag