


Anti-Ndufs2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - 075306

Item Number: CM0014766



Introduction

KO validated rabbit polyclonal antibody specific to NDUFS2. Reacts with human, mouse, rat. Suitable for WB, IHC-P, IF/ICC, ELISA applications. Ideal for mitochondrial complex I research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NDUFS2 Rabbit pAb
Remarks/Alias	CI-49; MC1DN6; S2
Species	Human
GeneID (human)	4720
GeneID	4720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 204-463 of human NDUFS2 (NP_004541.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075306
Protein Weight	53kDa
Shipping	Ice bag