


Anti-Ndufs3 Monoclonal Antibody For Wb, Ihc-P, Elisa - O75489

Item Number: CM0012257



Introduction

Rabbit monoclonal anti-NDUF53 antibody, validated for WB, IHC-P, ELISA. Recognizes human and rat NDUF53, a core subunit of mitochondrial Complex I (~30 kDa). Ideal for mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	NDUF53 Rabbit mAb
Remarks/Alias	CI-30; MC1DN8; NDUF53
Species	Human
Gene ID (Human)	4722
Gene ID	4722
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 165-264 of human NDUF53 (O75489).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O75489
Protein Weight	30 kDa
Shipping	Ice bag