


Anti-Ndufs3 Polyclonal Antibody For Wb And Elisa - O75489

Item Number: CM0026496



Introduction

Rabbit polyclonal antibody against human NDUF53, a core subunit of mitochondrial Complex I. Validated for WB and ELISA, with cross-reactivity to human and mouse. Ideal for studying mitochondrial respiratory chain function.

[Learn More](#)

Parameter	Value
Product Name	NDUF53 Rabbit pAb
Remarks/Alias	CI-30; MC1DN8; NDUF53
Species	Human
GeneID (Human)	4722
GeneID	4722
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-264 of human NDUF53 (NP_004542.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75489
Protein Weight	30kDa
Shipping	Ice bag