

Anti-Ndufs3 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - O75489

Item Number: CM0027432



Introduction

NDUF53 Rabbit Polyclonal Antibody validated for WB, IP, ELISA in human, mouse, rat. Specific to mitochondrial complex I subunit, ideal for respiratory chain research and diagnostics.

[Learn More](#)

Parameter	Value
Product Name	NDUF53 Rabbit pAb
Remarks/Alias	CI-30; MC1DN8; NDUF53
Species	Human
Gene ID (Human)	4722
Gene ID	4722
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-264 of human NDUF53 (NP_004542.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75489
Protein Weight	30kDa
Shipping	Ice bag