


Anti-Ndufs4 Monoclonal Antibody For Wb - O43181

Item Number: CM0011670



Introduction

Rabbit monoclonal antibody against human NDUF54 validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for mitochondrial complex I and respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUF54 Rabbit mAb
Remarks/Alias	AQDQ; CI-18; MC1DN1; CI-AQDQ; CI-18 kDa; NDUF54
Species	Human
Gene ID	4724
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human NDUF54 (O43181)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43181
Protein Weight	20kDa
Shipping	Ice bag