


Anti-Ndufs4 Rabbit Monoclonal Antibody For Wb, Elisa - O43181

Item Number: CM0007655



Introduction

Rabbit monoclonal antibody targeting human NDUF54 (Complex I-18 kDa), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial dysfunction and Complex I research.

[Learn More](#)

Parameter	Value
Product Name	NDUF54 Rabbit mAb
Remarks/Alias	AQDQ; CI-18; MC1DN1; CI-AQDQ; CI-18 kDa; NDUF54
Species	Human
Gene ID (Human)	4724
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-175 of human NDUF54 (NP_002486.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43181
Protein Weight	20kDa
Shipping	Ice bag