


Anti-Ndufs4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43181

Item Number: CM0016194



Introduction

Rabbit polyclonal antibody targeting human NDUF54, a mitochondrial Complex I subunit. Validated for WB, IHC-P, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Ideal for mitochondrial disease and respiratory chain studies.

[Learn More](#)

Parameter	Value
Product Name	NDUF54 Rabbit pAb
Remarks/Alias	AQDQ; CI-18; MC1DN1; CI-AQDQ; CI-18 kDa; NDUF54
Species	Human
GeneID (human)	4724
GeneID	4724
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human NDUF54 (NP_002486.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43181
Protein Weight	20kDa
Shipping	Ice bag