


Anti-Ndufs5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43920

Item Number: CM0024388



Introduction

Rabbit polyclonal antibody against human NDUF55 (CI-15kDa), a Complex I accessory subunit. Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for mitochondrial respiratory chain studies.

[Learn More](#)

Parameter	Value
Product Name	NDUF55 Rabbit pAb
Remarks/Alias	CI15K; CI-15k; NDUF55
Species	Human
Gene ID	4725 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human NDUF55 (NP_004543.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43920
Protein Weight	13kDa
Shipping	Ice bag