

Anti-Ndufs5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43920

Item Number: CM0020588



Introduction

Rabbit polyclonal antibody targeting human NDUFS5 (CI-15k), a Complex I accessory subunit. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Suitable for mitochondrial respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] NDUFS5 Rabbit pAb
Remarks/Alias	CI15K; CI-15k; [KD Validated] NDUFS5
Species	Human
Gene ID (Human)	4725
Gene ID	4725
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human NDUFS5 (NP_004543.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43920
Protein Weight	13kDa
Shipping	Ice bag