


Anti-Ndufs7 Rabbit Polyclonal Antibody For Wb, Elisa - 075251

Item Number: CM0018871



Introduction

Rabbit polyclonal antibody targeting human NDUFS7 (Complex I-20kD, PSST) for WB and ELISA. Essential for mitochondrial respiration studies. Validated in human samples; suitable for mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	NDUFS7 Rabbit pAb
Remarks/Alias	NDUFS7; CI-20; CI-20KD; MY017; PSST; NADH:ubiquinone oxidoreductase core subunit S7
Species	human
GeneID	374291
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 80-160 of human NDUFS7 (NP_077718.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	075251
Protein Weight	22kDa/23kDa
Shipping	Ice bag