


Anti-Ndufs8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O00217

Item Number: CM0013896



Introduction

Rabbit polyclonal antibody against NDUF58, a mitochondrial Complex I subunit, validated for WB, IF/ICC, ELISA in human and mouse. Ideal for studying mitochondrial respiratory chain function.

[Learn More](#)

Parameter	Value
Product Name	NDUF58 Rabbit pAb
Remarks/Alias	TYKY, CI-23k, CI23KD, MC1DN2, NDUF58
Species	Human
GeneID (Human)	4728
GeneID	4728
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 35-210 of human NDUF58 (NP_002487.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	O00217
Protein Weight	24kDa
Shipping	Ice bag