


Anti-Ndufv1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P49821

Item Number: CM0006899



Introduction

NDUFV1 rabbit monoclonal antibody (clone not specified) detects the 51 kDa mitochondrial complex I subunit across human, mouse, and rat samples. Validated for WB, IHC-P, and ELISA. Ideal for studies on electron transport chain, mitochondrial disease, and metabolomics.

[Learn More](#)

Parameter	Value
Product Name	NDUFV1 Rabbit mAb
Remarks/Alias	UQOR1; CI-51K; CI51KD; MC1DN4; NDUFV1
Species	Human
GeneID (Human)	4723
GeneID	4723
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 72-420 of human NDUFV1 (NP_009034.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49821
Protein Weight	51kDa
Shipping	Ice bag