


Anti-Ndufv1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49821

Item Number: CM0018031



Introduction

This rabbit polyclonal antibody against human NDUFV1 (NADH dehydrogenase flavoprotein 1, Complex I-51kD) is validated for western blot, IHC-paraffin, immunofluorescence/ICC, and ELISA; it cross-reacts with mouse and rat NDUFV1.

[Learn More](#)

Parameter	Value
Product Name	NDUFV1 Rabbit pAb
Remarks/Alias	UQOR1; CI-51K; CI51KD; MC1DN4; NDUFV1
Species	Human
GeneID (Human)	4723
GeneID	4723
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P49821
Protein Weight	51kDa
Shipping	Ice bag