


Anti-Ndufv1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49821

Item Number: CM0029635



Introduction

Rabbit polyclonal antibody targeting NDUFV1, a 51kDa core subunit of mitochondrial Complex I. Validated for WB, IHC-P, and ELISA with reactivity across human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	NDUFV1 Rabbit pAb
Remarks/Alias	UQOR1; CI-51K; CI51KD; MC1DN4; NDUFV1
Species	Human
GeneID (Human)	4723
GeneID	4723
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 72-420 of human NDUFV1 (NP_009034.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49821
Protein Weight	51kDa
Shipping	Ice bag