


Anti-Ndufv2 Polyclonal Antibody For Wb, Elisa - P19404

Item Number: CM0030327



Introduction

Rabbit polyclonal antibody against human NDUFV2 (CI-24k). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. NDUFV2 is a core subunit of mitochondrial respiratory chain Complex I, involved in electron transfer and implicated in mitochondrial disease. - P19404

[Learn More](#)

Parameter	Value
Product Name	NDUFV2 Rabbit pAb
Remarks/Alias	CI-24k; MC1DN7; NDUFV2
Species	Human
Gene ID (Human)	4729
Gene ID	4729
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-249 of human NDUFV2 (NP_066552.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19404
Protein Weight	27kDa
Shipping	Ice bag