

Anti-Necab1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N987

Item Number: CM0013089

Introduction

NECAB1 (EFCBP1) rabbit polyclonal antibody raised against amino acids 1-351. Suitable for western blot and ELISA, with confirmed cross-reactivity to human, mouse, and rat. Ideal for studying neuronal calcium-binding protein 1 in brain tissue.

[Learn More](#)

Parameter	Value
Product Name	NECAB1 Rabbit pAb
Remarks/Alias	EFCBP1; STIP-1; NECAB1
Species	Human
Gene ID	64168
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human NECAB1 (NP_071746.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N987
Protein Weight	41kDa
Shipping	Ice bag