

Anti-Necab3 Rabbit Polyclonal Antibody For Wb And Elisa - Q96P71

Item Number: CM0016811

Introduction

Rabbit polyclonal antibody targeting human NECAB3, suitable for WB and ELISA, with cross-reactivity to mouse and rat. A valuable tool for investigating amyloid-beta formation and Alzheimer's disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	NECAB3 Rabbit pAb
Remarks/Alias	NIP1; XB51; STIP3; EFCBP3; SYTIP2; APBA2BP; dj63M2.4; dj63M2.5; NECAB3
Species	Human
Gene ID	63941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-362 of human NECAB3 (NP_112508.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96P71
Protein Weight	44kDa
Shipping	Ice bag