

Anti-Necap2 Polyclonal Antibody For Wb, Elisa - Q9Nvz3

Item Number: CM0018122

Introduction

Rabbit polyclonal antibody targeting human NECAP2 (Adaptin ear-binding coat-associated protein 2). Validated for WB and ELISA, cross-reacts with mouse. Suitable for endocytosis and vesicle trafficking research. - Q9NVZ3.

[Learn More](#)

Parameter	Value
Product Name	NECAP2 Rabbit pAb
Remarks/Alias	NECAP2
Species	Human
GeneID (Human)	55707
GeneID	55707
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 194-263 of human NECAP2 (NP_060560.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NVZ3
Protein Weight	28kDa
Shipping	Ice bag