

Anti-Nectin 2/Cd112 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92692

Item Number: CM0011461



Introduction

Rabbit monoclonal antibody targeting human Nectin-2 (CD112, PVRL2) for WB, IHC-P, ELISA. Detects human, mouse, and rat protein; involved in T-cell modulation and viral entry.

[Learn More](#)

Parameter	Value
Product Name	Nectin 2/CD112 Rabbit mAb
Remarks/Alias	HVEB; PRR2; CD112; PVRL2; PVRR2; Nectin 2/CD112
Species	Human
GeneID (human)	5819
GeneID	5819
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92692
Protein Weight	58kDa
Shipping	Ice bag