

Anti-Nectin 2/Cd112 Rabbit Polyclonal Antibody For Wb And Elisa - Q92692

Item Number: CM0029191



Introduction

High-quality rabbit polyclonal antibody against human Nectin-2 (CD112, HVEB), validated for WB and ELISA. Ideal for T-cell signaling, cell adhesion, and viral entry research.

[Learn More](#)

Parameter	Description
Product Name	Nectin 2/CD112 Rabbit pAb
Remarks/Alias	HVEB; PRR2; CD112; PVRL2; PVRR2; Nectin 2/CD112
Species	Human
GeneID (Human)	5819
GeneID	5819
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 382-479 of human Nectin 2/CD112 (NP_002847.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92692
Protein Weight	58kDa
Shipping	Ice bag