

Anti-Nectin-4/Pvrl4 Rabbit Monoclonal Antibody For Fc, Elisa - Q96Ny8

Item Number: CM0009137



Introduction

Recombinant rabbit monoclonal antibody targeting human Nectin-4/PVRL4. Validated for FC and ELISA. Suitable for cancer and viral infection research, specifically measles virus entry studies.

[Learn More](#)

Parameter	Value
Product Name	Nectin-4/PVRL4 Rabbit mAb
Remarks/Alias	LNIR; PRR4; EDSS1; PVRL4; nectin-4
Species	Human
Gene ID	81607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-146 of human Nectin-4/PVRL4 (NP_112178.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	Q96NY8
Protein Weight	24kDa/55kDa
Shipping	Ice bag