


Anti-Nedd1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nhv4

Item Number: CM0030106



Introduction

Rabbit polyclonal antibody targeting human NEDD1, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Recognizes a 72kDa band. Ideal for mitosis progression and microtubule nucleation studies.

[Learn More](#)

Parameter	Value
Product Name	NEDD1 Rabbit pAb
Remarks/Alias	GCP-WD; TUBGCP7; NEDD1
Species	Human
GeneID (Human)	121441
GeneID	121441
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 411-660 of human NEDD1 (NP_690869.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NHV4
Protein Weight	72kDa
Shipping	Ice bag