

Anti-Nedd4 Monoclonal Antibody For Wb, Elisa - P46934

Item Number: CM0010266



Introduction

Anti-NEDD4 rabbit monoclonal antibody validated for WB and ELISA. Targets human NEDD4 (149kDa E3 ubiquitin ligase), relevant for ubiquitination, autophagy, and viral infection studies.

[Learn More](#)

Parameter	Value
Product Name	NEDD4 Rabbit mAb
Remarks/Alias	RPF1; NEDD4-1; NEDD4
Species	Human
GeneID (Human)	4734
GeneID	4734
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P46934
Protein Weight	149kDa
Shipping	Ice bag