

Anti-Nedd4 Polyclonal Antibody For Wb, Ihc-P, Elisa - P46934

Item Number: CM0015562



Introduction

Rabbit polyclonal antibody against human NEDD4, an E3 ubiquitin-protein ligase involved in receptor degradation and autophagy. Validated for WB, IHC-P, ELISA. Ideal for ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	NEDD4 Rabbit pAb
Remarks/Alias	RPF1; NEDD4-1; NEDD4
Species	Human
GeneID (Human)	4734
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-460 of human NEDD4 (NP_940682.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P46934
Protein Weight	149kDa
Shipping	Ice bag