

Anti-Nedd4L Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Pu5

Item Number: CM0027668



Introduction

NEDD4L Rabbit Polyclonal Antibody (Q96PU5) validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, and rat NEDD4L, an E3 ubiquitin ligase regulating TGF- β , autophagy, and ion channels. Ideal for signaling and ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	NEDD4L Rabbit pAb
Remarks/Alias	RSP5; PVNH7; NEDD4-2; NEDD4.2; hNEDD4-2; NEDD4L
Species	Human
GeneID (Human)	23327
GeneID	23327
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 836-955 of human NEDD4L (NP_056092.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PU5
Protein Weight	112kDa
Shipping	Ice bag