

Anti-Nedd4L Rabbit Monoclonal Antibody For Wb, Elisa - Q96Pu5

Item Number: CM0012441



Introduction

NEDD4L Rabbit mAb for WB and ELISA. Validated in human samples, this monoclonal antibody detects E3 ubiquitin ligase NEDD4L (~112 kDa). Sourced for research on ubiquitin pathways, channel regulation, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	NEDD4L Rabbit mAb
Remarks/Alias	RSP5; PVNH7; NEDD4-2; NEDD4.2; hNEDD4-2; NEDD4L
Species	Human
Gene ID (Human)	23327
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 876-975 of human NEDD4L (Q96PU5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96PU5
Protein Weight	112kDa
Shipping	Ice bag