

Anti-Nedd8 Rabbit Monoclonal Antibody For Wb, Elisa - Q15843

Item Number: CM0007516



Introduction

NEDD8 Rabbit Monoclonal Antibody validated for Western blot and ELISA, suitable for human, mouse, and rat samples. Targets ubiquitin-like protein NEDD8 involved in cell cycle control and protein degradation.

[Learn More](#)

Parameter	Value
Product Name	NEDD8 Rabbit mAb
Remarks/Alias	NEDD-8; NEDD8
Species	Human
Gene ID (Human)	4738
Gene ID	4738
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15843
Protein Weight	9kDa
Shipping	Ice bag