

Anti-Nedd8 Rabbit Polyclonal Antibody For Wb And Elisa - Q15843

Item Number: CM0017111



Introduction

Rabbit polyclonal antibody targeting human NEDD8, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Suitable for ubiquitin-like protein and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	NEDD8 Rabbit pAb
Remarks/Alias	NEDD-8; NEDD8
Species	Human
GeneID (Human)	4738
GeneID	4738
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-81 of human NEDD8 (NP_006147.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15843
Protein Weight	9kDa
Shipping	Ice bag