


Anti-Nedd9 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14511

Item Number: CM0022240



Introduction

Rabbit polyclonal antibody against human NEDD9 (HEF1/CAS-L). Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for cell adhesion, migration, and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	NEDD9 Rabbit pAb
Remarks/Alias	CAS2; CASL; HEF1; CAS-L; CASS2; NEDD9
Species	Human
GeneID (Human)	4739
GeneID	4739
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 390-560 of human NEDD9 (NP_006394.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14511
Protein Weight	93kDa
Shipping	Ice bag