


Anti-Nek2 Monoclonal Antibody For Wb, Elisa - P51955

Item Number: CM0012734



Introduction

Rabbit monoclonal antibody against human NEK2, validated for Western blot and ELISA. Detects human, mouse, and rat NEK2 (52kDa). Ideal for cell cycle and centrosome research.

[Learn More](#)

Parameter	Value
Product Name	NEK2 Rabbit mAb
Remarks/Alias	NLK1; RP67; NEK2A; HsPK21; PPP1R111; NEK2
Species	Human
GeneID (Human)	4751
GeneID	4751
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NEK2 (P51955)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51955
Protein Weight	52kDa
Shipping	Ice bag