

Anti-Nek3 Polyclonal Antibody For Wb And Elisa - P51956

Item Number: CM0029646

Introduction

Rabbit polyclonal antibody against human NEK3 kinase (P51956) validated for Western blot and ELISA, with reactivity in human, mouse, and rat. Suitable for neuronal morphogenesis and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	NEK3 Rabbit pAb
Remarks/Alias	HSPK36; NEK3
Species	Human
GeneID (Human)	4752
GeneID	4752
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 267-506 of human NEK3 (NP_689933.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51956
Protein Weight	58kDa
Shipping	Ice bag