

Anti-Nek6 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Hc98

Item Number: CM0009672

Introduction

This NEK6 Rabbit Monoclonal Antibody detects endogenous NEK6 in human, mouse, and rat samples by WB, IHC-P, ELISA. It targets a core kinase domain region and supports cell cycle, mitosis, and DNA damage research.

[Learn More](#)

Parameter	Value
Product Name	NEK6 Rabbit mAb
Remarks/Alias	SID6-1512; NEK6
Species	Human
Gene ID	10783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human NEK6 (Q9HC98)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC98
Protein Weight	36kDa
Shipping	Ice bag