


Anti-Nek7 Polyclonal Antibody For Wb, Elisa - Q8Tdx7

Item Number: CM0021233



Introduction

Rabbit polyclonal antibody against human NEK7, suitable for WB and ELISA. Reacts with human, mouse, and rat. Useful for studying mitotic cell cycle progression and NLRP3 inflammasome activation.

[Learn More](#)

Parameter	Value
Product Name	NEK7 Rabbit pAb
Species	Human
Gene ID (Human)	140609
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 257-301 of human NEK7 (NP_598001.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDX7
Protein Weight	35kDa
Shipping	Ice bag