


Anti-Nek7 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Tdx7

Item Number: CM0002629



Introduction

NEK7 Rabbit Monoclonal Antibody validated for WB and ELISA; reacts with Human, Mouse, Rat; targets Serine/threonine-protein kinase Nek7 involved in mitosis and NLRP3 inflammasome activation. Ideal for cell cycle and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	NEK7 Rabbit mAb
Remarks/Alias	NIMA-related kinase 7; NEK7
Species	Human
Gene ID (Human)	140609
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 203-302 of human NEK7 (Q8TDX7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDX7
Protein Weight	35kDa
Shipping	Ice bag