

Anti-Nek8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q86Sg6

Item Number: CM0014532

Introduction

High-quality rabbit polyclonal NEK8 antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat NEK8, a serine/threonine kinase critical for ciliary function and renal tubular integrity. Ideal for nephrology and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	NEK8 Rabbit pAb
Remarks/Alias	JCK, NPHP9, RHPD2, NEK12A, NEK8
Species	Human
GeneID (Human)	284086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-390 of human NEK8 (NP_835464.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86SG6
Protein Weight	75kDa
Shipping	Ice bag