

Anti-Nek9 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Td19

Item Number: CM0021144



Introduction

NEK9 rabbit polyclonal antibody for WB, IHC-P, and ELISA applications. Detects human, mouse, and rat NEK9 (~107 kDa). Ideal for studying mitotic progression and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	NEK9 Rabbit pAb
Remarks/Alias	NC; APUG; NERCC; LCCS10; NERCC1; NEK9
Species	Human
Gene ID (Human)	91754
Gene ID	91754
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 730-979 of human NEK9 (NP_149107.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TD19
Protein Weight	107kDa
Shipping	Ice bag