


## Anti-Nelfcd Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Ixx7

Item Number: CM0025107



### Introduction

Anti-NELFCD polyclonal antibody, reacting with human, mouse, rat NELFCD (TH1/TH1L). Validated for WB, IHC-P, ELISA. Suitable for transcription elongation and NELF complex studies.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | NELFCD Rabbit pAb                                                                                                  |
| Remarks/Alias    | TH1; TH1L; NELF-C; NELF-D; HSPC130; NELFCD                                                                         |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 51497                                                                                                              |
| GeneID           | 51497                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 40-182 of human NELFCD (NP_945327.3) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q8IXH7                                                                                                             |
| Protein Weight   | 66kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |