


Anti-Nelfe Polyclonal Antibody For Wb, Ihc-P, Elisa - P18615

Item Number: CM0028798



Introduction

Rabbit polyclonal antibody against human NELFE (RD, RDBP), validated in WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. NELFE is a transcription elongation repressor and component of the NELF complex.

[Learn More](#)

Parameter	Value
Product Name	NELFE Rabbit pAb
Remarks/Alias	RD; RDP; RDBP; D6S45; NELF-E; NELFE
Species	Human
Gene ID (Human)	7936
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NELFE (NP_002895.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18615
Protein Weight	43kDa
Shipping	Ice bag