

Anti-Nell2 Polyclonal Antibody For Wb, Elisa - Q99435

Item Number: CM0019632



Introduction

Rabbit polyclonal antibody against human NELL2 (Protein kinase C-binding protein NELL2). Validated for WB and ELISA applications. Cross-reacts with mouse and rat. Ideal for neural tissue and axon guidance research.

[Learn More](#)

Specification	Detail
Product Name	NELL2 Rabbit pAb
Remarks/Alias	NRP2; NELL2
Species	Human
Gene ID (Human)	4753
Gene ID	4753
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-300 of human NELL2 (NP_006150.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99435
Protein Weight	91kDa
Shipping	Ice bag