

Anti-Nephrin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Qzs7

Item Number: CM0004146

Introduction

High-specificity Nephrin Rabbit Monoclonal Antibody validated for WB, IF-P, IHC-P, and ELISA in mouse. Targets podocyte slit diaphragm protein; essential for kidney filtration and nephropathy research. UniProt: Q9QZS7.

[Learn More](#)

Parameter	Value
Product Name	Nephrin Rabbit mAb
Remarks/Alias	nephrin; NephrinB
Species	Mouse
GenelD (human)	4868
GenelD	54631
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q9QZS7
Protein Weight	135 kDa/136 kDa
Shipping	Ice bag